




2900 12<sup>th</sup> Avenue North, Suite 3E  
Lower Level Yellowstone Medical Center East  
Billings, Montana 59101  
Tel. 406.237.5525  
Fax. 406.237.5530  
www.bigskydiagnosticimaging.com

## **Viewing Patient Images on the Fuji Synapse from Big Sky Diagnostic Imaging of Billings**


1. Double Click on the Synapse icon on the desktop



2. Type in your user name and password. (the user name will have to be proceeded by “bsdi\”, for example, if the user name is “jsmith” then the user name when logging in would be “bsdi\jsmith”)
3. Double click on the “Big Sky Diagnostic Imaging” name
4. Double click on the “All Studies” folder listed



5. Type the complete last name of the patient under the empty box below “Patient name”



Synapse/Mercury/All Studies					
Patient Name	Accession No	Patient ID	Proc Description	Study Date Time	Mod...
b					
Bucknell, Steve	D980076	FMSU0076	System Test St...	11/19/1998 9:3...	CR
Bucknell, Steve	D980076-1	FMSU0076	System Test St...	11/27/2002 9:4...	CR
Busse, Ann C	2382294	000221712	CT ABDOMEN	5/30/2000 11:2...	CT
Bbbbbbbbbbbbbb, ...	1234567890...	AAAAAAAAAA...	UNKNOWN	9/6/2001 6:07:...	CT
Bell, Randy U	1213164711...	285590055	Exercise Echoc...	12/14/2002 4:4...	US

6. When patient, and correct study appear, double click on patient name

Synapse/Thender/Unread Studies			
Patient Name	Accession No	Patient ID	Proc Description
Duncan,, Michael	D99062	FMSU021	ABDOMEN, AB...
Deets,, Eric	D99073	FMSU033	PELVIS, HIP JO...
Zoller,, Katherina	D99078	FMSU038	ABDOMEN, UP...
Feingold, Steve	D980090	FMSU0085	ABDOMEN, AB...


7. The study will now open and be displayed (this may take a moment depending on your computers speed)

8. The images can be scrolled through by clicking on the desired sequence with one left mouse click. You can then use your middle scroll on your mouse or your arrow keys on your keyboard.

9. To view the report you can click on the patient power jacket icon



You then click on the “DOCS” icon within the power jacket. You will see the patient name appear within the documents listed. Double click on the patient name and report will open up.



The screenshot shows a software interface with a sidebar on the left containing icons for 'Home', 'Files', 'Notes', 'Print', and 'Print All'. The main area displays patient information: 'Patient Name: Michael Duncan', 'Accession #: D99062', 'Patient ID: FMSU021', and 'Proc Description: ABDOMEN, AB...'. Below this is a list of documents. A red arrow points to the bottom of this list, which contains several entries with document names and icons.

10. You can access the tools for viewing with one “right” click inside the image box. (i.e.: Windowing, etc.)
11. When you have finished viewing the case and/or are done utilizing Synapse, close out by clicking on the “X”
12. Right Click on the Synapse icon and “**log off**” of Synapse



**MONTANA'S ONLY ULTRA-HIGH-FIELD 3.0 TESLA MRI SYSTEM**  
*World Class Technology / Friendly Outpatient Setting*